



European Directorate for the Quality of Medicines & HealthCare
European Pharmacopoeia (Ph. Eur.)
7, Allée Kastner CS 30026, F-67081 Strasbourg (France)
Tel. +33 (0)3 88 41 20 35 Fax. + 33 (0)3 88 41 27 71
For any questions: www.edgm.eu (HelpDesk)

Ph. Eur. Reference Standard - LEAFLET

BISOPROLOL FOR SYSTEM SUITABILITY METHOD A CRS batch 1

Intended use

This leaflet supplements the currently valid European Pharmacopoeia monograph(s) and/or general chapter(s) describing the suitable use of this Reference Standard.

Further information about the Reference Standards is available in the on-line catalogue currently at http://crs.edqm.eu (such as, batch validity statements and material safety data sheets).

Instruction

Bisoprolol for system suitability method A CRS 1 has been renamed as bisoprolol for peak identification CRS 1.


A chromatogram obtained with bisoprolol for system suitability method A CRS 1 in accordance with the monograph 04/2008:1710 is provided in the annex 1.

A chromatogram obtained with bisoprolol for peak identification CRS 1 in accordance with the monograph 01/2012:1710 is provided in the annex 2.

* 0 0 Y 0 0 0 0 8 7 7 0 4 . 0 *

Catalogue code: Y0000877 Revision 04 Date of issue: 20/05/2011 1/3

EUROPEAN DIRECTORATE FOR THE QUALITY OF MEDICINES



Project Name: LC05165BisoproloISST Result Id 1591




Catalogue code: Y0000877 Revision 04 Date of issue: 20/05/2011 2/3

Annex 2



LIQUID CHROMATOGRAPHY REPORT

Bisoprolol for Peak Identification CRS1



Project Name: LC06792BISOPROLOL FUMARATE

Result ld 1605



Catalogue code: Y0000877 Revision 04 Date of issue: 20/05/2011 3/3